

## >> Multifunction videosystem



>> The consolidated unit  
for the maximal user friendliness




## >> Multifunction videosystem

### Product profile:

- consolidated system consisting of a light source, camera, TFT-monitor and (optional) archiving (USB)
- the easiest use
- interesting design
- high quality manufacture – „made in Germany“

### User-friendliness along with the maximal scope of operation:

- control functions: white balance, window, boost, Anti-Moiré
- can be combined with light guides of different lengths and diameters



### Technical specifications:

- video resolution 470 TV-lines
- shutter range 1/50 sec – 1/100.000 sec
- sensitivity 0,9 lx (f:1,4)
- video outputs 1 x Y/C (S-VHS), 1 x BNC
- lamps xenon metal halide short arc lamp 24 or 50 W, 5.500K colour temperature
- power supply 90 VAC – 260 VAC
- weight CCU: 6kg  
camera head: 65g
- class 1
- manufactured and tested acc. IEC 601-1 / CE  
EN 60601-1  
93 / 42 ECC

### Scope of delivery:

Multifunction videosystem + AC/DC adapter cable + S-VHS-cable + BNC-cable

### Optional accessories:

light adjustment device (stepless, mechanical)