


>> mobiLED⁺



>> The LED light source
with maximum power



Product profile:

- unrivalled illuminating power
- extraordinary compact design
- stepless brightness adjustment

User-friendliness along with the maximal scope of operation:

- rechargeable lithium ion battery enables battery-supplied operation
- ultimate compatibility assured by **ACMI connection**
- no exchange of lamps necessary due to LED technique
- 45 minutes non-stop operation at maximum power
- **waterproof**, completely submersible
- 1 year warranty



Technical specifications:

- | | |
|--|--|
| ■ light output | minimum 100 lm |
| ■ charge time | 2 h |
| ■ weight | 90 g |
| ■ battery type | lithium ion NP-700, 3.6V, 1,000mAh minimum |
| ■ safe operating ambient temperature range | +15°C – +30°C |
| ■ safe storage and transport temperature range | -25°C – +50°C |
| ■ relative humidity | 0 % – 95 % |

Scope of delivery:

mobiLED⁺ + lithium ion battery + battery charger + AC/DC adapter cable

