

LED Dual Gooseneck Lights Illuminator with 2pcs 6 WATT LED Gooseneck Lights

Features:

Step-less dimming

Easy change detachable

Gooseneck light source

Eleven different colour options

High index of refraction glass lense



Specification:

Size: 140mm x 106mm x 650mm (L x W x L)

Neck Length: 590mm

Gooseneck Length: 450mm

Power: 6W (2 * 3W)

LED Colour option: White / Blue (440nm wavelength) / Green (530nm wavelength) / UV365 (365nm wavelength) / UV395 (395nm wavelength) / UV420 (420nm wavelength)

Input Voltage: 85V-265V

Custom made Light Guide Cable adapters available for Endoscopes



Contact: VisionScope Technologies Ltd. Tel: 0845 299 7440.

Email: visionscope@btinternet.com. Web: www.vst.org.uk