



# Kevin Sebley - Endoscopy



## LAPAROSCOPY



>> Laparoscopy.

Laparoscopes		l= 300 mm	l= 310 mm
Ø 5,0 mm	0°	30-0071-00	
Ø 5,0 mm	30°	30-0072-00	
Ø 6,5 mm	0°	30-0173-00	
Ø 6,5 mm	30°	30-0074-00	
Ø 10 mm	0°		30-0075-00
Ø 10 mm	30°		30-0076-00

autoclavable (134°C / 273°F)

Children Laparoscopes		l= 200 mm
Ø 4,0 mm	0°	30-0077-00
Ø 4,0 mm	30°	30-0078-00

autoclavable (134°C / 273°F)

Mini Laparoscopes		l= 300 mm
Ø 2,0 mm	0°	30-0079-00

scalable

Children Mini Laparoscopes		l= 200 mm
Ø 2,0 mm	0°	30-0136-00

scalable

Thorascopes		l= 170 mm
Ø 10 mm	0°	30-0082-00
Ø 10 mm	30°	30-0083-00

autoclavable (134°C / 273°F)

# Kevin Sebley - Endoscopy




**Web: [www.arthroscopyuk.co.uk](http://www.arthroscopyuk.co.uk)**

**Email: [visionscope@live.co.uk](mailto:visionscope@live.co.uk)**

**Tel: 07761 674499**

Hipp Endoskop Service GmbH



>> your requirements are our aims.