


HIPP Endoskop Service GmbH


>> your requirements are our aims.



Technische Produktinformation Hysteroskope
Technical Product Information Hysteroscopes



Hysteroskope Hysteroscopes



Hysteroskope – Storz-kompatibel
Hysteroscopes – Storz-compatible

autoklavierbar (134°C / 273°F)
autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
2.7	0°	302	x	30-0061-00
2.7	12°↑	302	x	30-0124-00
2.9	0°	302	x	30-0142-00
2.9	30°↑	302	x	30-0144-00
4.0	0°	302	x	30-0063-00
4.0	12°↑	302	x	30-0064-00
4.0	30°↑	302	x	30-0065-00
4.0	70°↑	302	x	30-0066-00

Kevin Sebley - Endoscopy

4 Beatrice Road, Wisbech, Cambridgeshire, PE13 3PE

Tel: 07761 674499

Email: visionscope@live.co.uk

Kevin Sebley, Product Consultant & Engineer

