


# HIPP Endoskop Service GmbH


>> *your requirements are our aims.*



Technische Produktinformation Lichtkabel  
Technical Product Information light cables



**Lichtkabel**  
light cables




- Kaltlichtkabel mit Durchmessern von 3,5 mm oder 4,8 mm
- Arbeitslängen von 1800 mm bis 3000 mm
- für Xenon und LED Lichtquellen geeignet
- mit innenliegendem Knickschutz und Glasstab für Temperaturbeständigkeit
- autoklavierbar 134°C/273°F
- Sonderlängen möglich
- cold light cable with diameter of 3.5 mm or 4.8 mm
- working length from 1800 mm to 2300 mm
- suitable for xenon and LED light sources
- with internal kink protection and glass rod for temperature resistance
- autoclavable 134°C/273°F
- Special lengths available

**Lichtkabel**  
light cables







Ø [mm]	Nutzlänge [mm]	autoklavierbar (134°C / 273°F)	Artikel-Nr.
Ø [mm]	working length [mm]	autoclavable (134°C / 273°F)	article-no.
3,5	1800	✓	14-0037-23
3,5	2300	✓	14-0038-23
3,5	3000	✓	14-0039-23
4,8	1800	✓	14-0043-23
4,8	2300	✓	14-0044-23
4,8	2500	✓	14-0045-23
4,8	3000	✓	14-0047-23

**Lichtleiter-Adapter-System  
Light guide adapter system**



**Adapter Endoskop-Seite  
Adapters optic side**

System system		Artikel-Nr. article-no.
Storz, Olympus, Dyonics, Linvatec, Atmos, Schöllly, Asap, Stryker, Aesculap		13-1034-00
Wolf "snap"		13-1036-00
Wolf "clip"		13-1050-00
ACMI, British		14-1065-00

System system		Artikel-Nr. article-no.
Storz		13-1031-00
Olympus		13-1032-00
Wolf		13-1033-00
ACMI/British, Stryker		14-0054-00
Olympus, ACM		14-0060-00
Schott		14-0062-00
Olympus special		14-0063-00

## Kevin Sebley - Endoscopy

4 Beatrice Road, Wisbech, Cambridgeshire, PE13 3PE

Tel: 07761 674499

Email: [visionscope@live.co.uk](mailto:visionscope@live.co.uk).

Kevin Sebley, Product Consultant & Engineer

