


HIPP Endoskop Service GmbH


>> *your requirements are our aims.*



Technische Produktinformation Sinuskope
Technical Product Information Sinusscopes



Sinuskope Sinusscopes



Sinuskope – Storz-kompatibel
Sinusscopes – Storz-compatible

autoklavierbar (134°C / 273°F)
autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
2.7	0°	110	x	30-0029-00
2.7	30°	110	x	30-0030-00
2.7	70°	110	x	30-0031-00
2.7	0°	175	x	30-0032-00
2.7	30°	175	x	30-0033-00
2.7	70°	175	x	30-0034-00
2.7	0°	187.5	x	30-0035-00
2.7	30°	187.5	x	30-0036-00
2.7	70°	187.5	x	30-0037-00
4.0	0°	175	x	30-0038-00
4.0	30°	175	x	30-0039-00
4.0	45°	175	x	30-0028-00
4.0	70°	175	x	30-0040-00

Kevin Sebley - Endoscopy

4 Bzatriez Road, Wisbech, Cambridgeshire, PE13 3PE

Tel: 07761 674499

Email: visionseope@livz.co.uk

Kevin Sebley, Product Consultant & Engineer

