



HIPP Endoskop Service GmbH

>> *your requirements are our aims.*



Technische Produktinformation Uretero-Renokope
Technical Product Information Uretero-Renoscope




Uretero-Renokope Uretero-Renoscopes




Marberger



Bichler



5° ↓



12° ↓

Uretero-Renokope Uretero-Renoscopes

autoklavierbar (134°C / 273°F)
autoclavable (134°C / 273°F)

Ø [Cr/Fr] Ø [Cr/Fr]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	Modell Type	Artikel-Nr. article-no.
4.5 / 6.5	5° ↓	315	Marberger	30-0511-00
4.5 / 6.5	5° ↓	315	Bichler	30-0513-00
4.5 / 6.5	5° ↓	430	Marberger	30-0510-00
4.5 / 6.5	5° ↓	430	Bichler	30-0512-00
6.0 / 7.5	12° ↓	315	Marberger	30-0558-00
6.0 / 7.5	12° ↓	315	Bichler	30-0552-00
6.0 / 7.5	5° ↓	315	Marberger	30-0517-00
6.0 / 7.5	5° ↓	315	Bichler	30-0519-00
6.0 / 7.5	12° ↓	430	Marberger	30-0514-00
6.0 / 7.5	12° ↓	430	Bichler	30-0528-00
6.0 / 7.5	5° ↓	430	Marberger	30-0518-00
6.0 / 7.5	5° ↓	430	Bichler	30-0520-00
8.0 / 9.8	12° ↓	315	Marberger	30-0069-00
8.0 / 9.8	12° ↓	315	Bichler	30-0395-00
8.0 / 9.8	12° ↓	430	Marberger	30-0070-00
8.0 / 9.8	12° ↓	430	Bichler	30-0403-00
8.5 / 11.5	12° ↓	315	Marberger	30-0404-00
8.5 / 11.5	12° ↓	315	Bichler	30-0406-00
8.5 / 11.5	12° ↓	430	Marberger	30-0405-00
8.5 / 11.5	12° ↓	430	Bichler	30-0407-00



