


# HIPP Endoskop Service GmbH


>> *your requirements are our aims.*



Technische Produktinformation Cystoskope  
Technical Product Information Cystoscopes




## Cystoskope Cystoscopes



Cystoskope – Storz-kompatibel  
Cystoscopes – Storz-compatible

autoklavierbar (134°C / 273°F)  
autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
2.7	0	187.5	x	30-0269-00
2.7	30	187.5	X	30-0268-00
2.7	0°	302	x	30-0055-00
2.7	30°	302	x	30-0057-00
2.7	70°	302	x	30-0054-00
2.9	0°	302	x	30-0145-00
2.9	30°	302	x	30-0147-00
4.0	0°	302	x	30-0055-00
4.0	12°	302	x	30-0056-00
4.0	30°	302	x	30-0057-00
4.0	70°	302	x	30-0058-00



Cystoskope – Olympus-kompatibel  
Cystoscopes – Olympus-compatible

autoklavierbar (134°C / 273°F)  
autoclavable (134°C / 273°F)

Ø [mm] Ø [mm]	Blickrichtung direction of view	Nutzlänge [mm] working length [mm]	HD	Artikel-Nr. article-no.
2.7	30°	282	x	30-0054-05

## Kevin Sebley - Endoscopy

4 Beatriex Road, Wisbech, Cambridgeshire, PE13 3PE

Tel: 07761 674499

Email: [visionscope@livz.co.uk](mailto:visionscope@livz.co.uk)

Kevin Sebley, Product Consultant & Engineer

